

Appendix B: PBN Sample Collection and Processing Recording Form

This template is intended for your convenience, only. This information must be submitted to PBN via the REDCap Sample Collection and Processing Form at <https://redcap.link/PBNSampleForm>. If you have any questions about sample processing for this study, visit psychiatricbiomarkers.org or email pbn@iu.edu.

PBN Subject ID: _____

Visit: BL 12M 24M 36M 48M 60M

Kit Number (6-7 digits): _____

Participant Type: SCZ BD CHR Control

Sex (genetic sex, not gender): M F Other

Phlebotomy Details:

Date of Blood Draw:
[MM/DD/YYYY]

Time of Draw:
24 HR clock

Fasting Status:

Blood Processing:

If processing was performed per SOP, checkbox may be used in place of text field, where provided.

10ml Serum (red-top) Tubes	Standard Procedures (See PBN MOP for details)
Total volume collected for serum (mL):	From one 10ml serum tubes; 10ml expected.
Time spin started (24 HR):	Tube should be incubated upright at room temperature for at least 30 (but <60) minutes prior to centrifugation.
Duration of centrifugation (min):	<input type="checkbox"/> SOP – spin for 15 minutes
Temperature of centrifuge (°C):	<input type="checkbox"/> SOP – spin at 4°C
Force of centrifugation (x g):	<input type="checkbox"/> SOP – spin at 1500 x g
Was Serum NMDA aliquots created?: <i>Use orange-capped cryovial</i>	One 1ml serum NMDA cryovial expected, <i>only if subject is SCZ, BD or CHR</i>
# of serum aliquots created: <i>Use red-capped cryovials</i>	Two-three 1.5ml serum cryovials expected. If low volume draw occurs, please generate as many 1.5ml aliquots as possible
Time aliquots frozen (24 HR):	Cryovials should be frozen upright within 2 hours of collection
Storage temperature (°C):	<input type="checkbox"/> SOP – store at -80°C (±10°) Samples <i>may</i> be stored at -20°C for 2-4 hours.

10ml NaHep (green-top) Tube		Standard Procedures (See PBN MOP for details)
Total volume collected for PBMCs (mL):		Tube collected at each visit. Tube should be shipped ambient on the day of collection

10ml EDTA (purple-top) tubes		Standard Procedures (See PBN MOP for details)
Total volume collected for plasma and buffy coat (mL):		From one 10ml EDTA tubes; 10ml expected
Time spin started (24 HR):		10ml EDTA tube should be spun within 30 minutes of collection
Duration of centrifugation (min):		<input type="checkbox"/> SOP – spin for 15 minutes
Temperature of centrifuge (°C):		<input type="checkbox"/> SOP – spin at 4°C
Force of centrifugation (x g):		<input type="checkbox"/> SOP – spin at 1500 x g
# of plasma aliquots created: <i>Use purple-capped cryovials</i>		From 10ml EDTA tube, three-four 1.5ml plasma cryovials expected. If low volume draw occurs, please generate as many 1.5ml aliquots as possible
# of buffy coat aliquots created: <i>Use clear-capped cryovials</i>		From 10ml EDTA tube, one buffy coat cryovial expected. Each 10ml EDTA tube will produce 1 buffy coat. Buffy coats are ~750ul
Time aliquots frozen (24 HR):		Cryovials should be frozen upright within 2 hours of collection
Storage temperature (°C):		<input type="checkbox"/> SOP – store at -80°C (±10°) Samples <i>may</i> be stored at -20°C for 2-4 hours.

2.5ml RNA PAXgene™ tubes		Standard Procedures (See PBN MOP for details)
Total volume collected for whole blood (mL):		From two 2.5ml RNA PAXgene™ tubes; 5ml expected
Number of tubes collected:		Two tubes expected
Time tubes frozen (24 HR):		Tubes should be frozen upright within 2 hours of collection
Storage temperature (°C):		<input type="checkbox"/> SOP – store at -80°C (±10°) Samples <i>may</i> be stored at -20°C for 2-4 hours.

Notes (including any deviations from protocol):

Lumbar Puncture Details:

Date of Blood Draw:

[MM/DD/YYYY]

Time of Draw:

24 HR clock

Fasting Status:

Duration of LP:

Minutes

CSF Processing:*If processing was performed per SOP, checkbox may be used in place of text field, where provided.*

Processing Details		Standard Procedures (See PBN MOP for details)
Total volume collected for CSF (mL):		15ml expected
Time spin started (24 HR):		CSF should be spun within 30 minutes of collection
Duration of centrifugation (min):		<input type="checkbox"/> SOP – spin for 10 minutes
Temperature of centrifuge (°C):		<input type="checkbox"/> SOP – spin at 4°C
Force of centrifugation (x g):		<input type="checkbox"/> SOP – spin at 300 x g
Was a CSF NMDA aliquot created? <i>Use orange-capped cryovial</i>		One 1ml CSF NMDA cryovial expected for BL visit, <i>only if subject is SCZ, BD or CHR</i>
# of CSF aliquots created: <i>Use clear-capped cryovials</i>		Twenty 0.5ml serum cryovials expected.
# of CSF Cell Pellet aliquots created: <i>Use orange-capped cryovial</i>		One 0.5ml CSF Cell Pellet cryovial expected
Time aliquots frozen (24 HR):		Cryovials should be frozen upright within 2 hours of collection
Storage temperature (°C):		<input type="checkbox"/> SOP – store at -80°C (±10°) Samples <i>may</i> be stored at -20°C for 2-4 hours.

Notes (including any deviations from protocol):

Specimen Collection and Processing Form

Please complete the Specimen Collection and Processing Form, below.

Study Site

- Icahn School of Medicine at Mount Sinai
- Perelman School of Medicine at UPenn
- Yale School of Medicine
- Northwell Health

Email address of staff member completing this form

Note: A copy of the completed sample form and the shipping manifests will be sent to this address.

PBN Subject ID:

Is this participant a case or control?

- SCZ
 - BD
 - CHR
 - Control
- (CHR allowed for follow-up ONLY)

Subject's biological sex (used for DNA quality control)

- Male
- Female
- Other

Visit

- Baseline Visit
- 12 Month Visit
- 24 Month Visit
- 36 Month Visit
- 48 Month Visit
- 60 Month Visit
- Unscheduled Visit

Kit Number

Blood Collection and Processing

Date of blood collection

Time of blood collection

(Use 24 Hour clock.)

Patient's fasting status at time of blood collection

- Fasted
 Followed low-fat diet
 Not fasted/no dietary limitation

1. SERUM (one red-top 10mL serum tube)

Was blood collected and processed for SERUM?

- Yes
 No

Blood volume collected for SERUM

(mL)

Reason volume was less than standard

- Difficult stick/poor veins
 Patient dehydrated
 Bad tube vacuum
 Other

Time of SERUM tube centrifugation

(Use 24 Hour clock.)

Rate of SERUM tube centrifugation

(x g)

Duration of SERUM tube centrifugation

(minutes)

Temperature of SERUM tube centrifugation

(degrees Celsius)

Total volume of SERUM collected

(mL)

Was a 1 mL aliquot of serum created for NMDA receptor testing?

- Yes
 No
(Required for all SCZ/BP/CHR visits.)

If yes, set this aliquot aside to be shipped to UPENN.

Number of SERUM aliquots created to ship to the repository at IU

(Aliquot to 1.5 mL, if possible.)

Was a residual SERUM aliquot (less than 1.5 mL) created?

- Yes
 No
(If YES, please ensure the residual aliquot is capped with a blue top.)
-

What is the approximate volume of the residual SERUM aliquot?

_____ (mL)

Time SERUM was placed in freezer

_____ (Use 24 Hour clock.)

SERUM storage temperature

_____ (degrees Celsius)

SERUM notes

2. PBMCs (two green-top 10mL sodium heparin tubes)

Confirm participant consented to PBMC collection:

- Yes
 No
-

Was blood collected for PBMCs?

- Yes
 No
(PBMCs should be collected at all visits.)
-

Number of tubes collected for PBMCs

Total blood volume collected for PBMCs

_____ (mL)

Reason volume was less than standard

- Difficult stick/poor veins
 Patient dehydrated
 Bad tube vacuum
 Other
-

PBMC notes

3. PLASMA (one purple-top 10mL EDTA tube)

Was blood collected and processed for PLASMA?

- Yes
 No
-

Blood volume collected for PLASMA

_____ (mL)

Reason volume was less than standard

- Difficult stick/poor veins
 Patient dehydrated
 Bad tube vacuum
 Other

Time of PLASMA tube centrifugation

(Use 24 Hour clock.)

Rate of PLASMA tube centrifugation

(x g)

Duration of PLASMA tube centrifugation

(minutes)

Temperature of PLASMA tube centrifugation

(degrees Celsius)

Total volume of PLASMA collected

(mL)

Number of PLASMA aliquots created

(Aliquot to 1.5 mL, if possible.)

Was a residual PLASMA aliquot (less than 1.5 mL) created?

- Yes
 No
(If YES, please ensure the residual aliquot is capped with a blue top.)

What is the approximate volume of the residual PLASMA aliquot?

(mL)

Was the BUFFY COAT collected?

- Yes
 No

Time PLASMA and BUFFY COAT were placed in freezer

(Use 24 Hour clock.)

PLASMA and BUFFY COAT storage temperature

(degrees Celsius)

PLASMA notes

4. RNA (two 2.5mL PAXGene™ tubes)

Was blood collected for RNA?

- Yes
 No

Number of PAXGene™ tubes collected for RNA

Blood volume collected for RNA

Note: Max collection volume per tube is 2.5 ml.

_____ (mL)

Reason volume was less than standard

- Difficult stick/poor veins
- Patient dehydrated
- Bad tube vacuum
- Other

Date RNA were frozen

Time RNA were placed in freezer

_____ (Use 24 Hour clock.)

RNA storage temperature

_____ (degrees Celsius)

RNA notes

CSF Collection and Processing

Was CSF collected? Yes
 No

Why was CSF not collected?

Date of CSF collection

Time of CSF collection

(Use 24 Hour clock.)

Approximate duration of CSF collection (in minutes)

Note: We are interested in the duration of the CSF collection, specifically the length of time the collection needle is in the subarachnoid space.

(minutes)

Patient's fasting status at time of CSF collection Fasted
 Followed low-fat diet
 Not fasted/no dietary limitation

Was CSF submitted for clinical labs? Yes
 No

Was CSF collected using a syringe(s)? Yes
 No

Note: Drip method is preferred.

Total volume of CSF collected (incl. CSF submitted for clinical labs)

(mL)

Reason volume collected was less than standard Difficult collection/patient physiology
 Patient dehydrated
 Other

BASELINE SCZ/BP/CHR ONLY: Was a 1 mL aliquot of CSF created for NMDA receptor testing? Yes
 No

If yes, set this aliquot aside to be shipped to UPENN.

Time of CSF tube centrifugation

(Use 24 Hour clock.)

Rate of CSF tube centrifugation

(x g)

Duration of CSF tube centrifugation

(minutes)

Temperature of CSF tube centrifugation

(degrees Celsius)

Number of CSF aliquots created to ship to repository
at IU

(Aliquot to 500uL, if possible.)

Was the CSF PELLETT collected?

Yes
 No

Was the CSF PELLETT frozen in a prepared Mr. Frosty™?

Yes
 No

Time CSF, CSF NMDA (if collected), and CSF PELLETT were
placed in freezer

(Use 24 Hour clock.)

CSF, CSF NMDA (if collected), and CSF PELLETT initial
storage temperature

(degrees Celsius)

Was the CSF PELLETT transferred to liquid nitrogen
storage?

Yes
 No
(Required if samples held on site >48 hours.)

CSF notes
